



Recombinant *Oryza sativa* subsp. *japonica*

Putative 14-3-3-like protein GF14-H (GF14H)

Product Code	CSB-MP650087OFG
Abbreviation	GF14H
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2R1D5
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MKEREKVVRL AKLAEQAERY DDMVEFMKTL ARMDVDMSAE ERLIFSVMGFK KTIGARRASW RILESLEQKV TAGDQPGVTI NGYKKKVEDE LRAVCNEVLS IIAHCLPLA NSGENVVFFY KMKGDYYRYL AEFSTGTEKK AATDQSLMAY QAWPCAQLFS LLEIMNSPER ASQVAKQALD EATAEINSAG VEGYKDSMLM MQLLENLAL WTSELTGGET SKDDDVMEG
Source	Mammalian cell
Target Names	GF14H
Protein Names	Recommended name: Putative 14-3-3-like protein GF14-H Alternative name(s): G-box factor 14-3-3 homolog H
Expression Region	1-230
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.