



Recombinant Human Putative macrophage-stimulating protein MSTP9 (MST1P9)

Product Code	CSB-MP650360HU
Abbreviation	MST1P9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2TV78
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	QPGQRSPSND FQVLRGTELQ HLLHAVVPGP WQEDVADAAEE CAGRCGPLMD CWAFHYNVSS HGCQLLPWTQ HSPHSRLWHS GRCDLFQEKG EWGYMPTLRN GLEENFCRNP DGDPGGPWCH TTDPAVRFQS CSIKSCRVA CVWCNGEEYR GAVDRTESGR ECQRWDLQHP HQHPFEPGKF LDQGLDDNYC RNPDGSERPW CYTTDPQIER EFCDLPRCGS EAQPRQEATS VSCFRGKGEG YRGTTANTTTA AYLASVGTRK SHISTDLRQK NTRASEVGGG AGVGTCCCGD LRENFCWNLD GSEAPWCFTL RPTGTRVGFCY QIRRCTDDVR PQDCYHGAGE QYRGTVSKTR KGVQCQRWSA ETPHKLQALT LGRHALMSGT RAWKWLRLPC HDFAPAPASV HIYLRTACTT GGELLDPDPDG DSHGPWCYTM DPRTPFDYCA LRRCDQVQFE KCGKRVDRLD QRRSKLRVAG GHPGNSPWTV SLRNRQGQHF CAGSLVKEQW ILTARQCFSS CHMPLTGYEV WLGTLFQNPQ HGEPGLQRVP VAKMLCGPSG SQLVLLKLER SVTLNQRVAL ICLPPEWYVV PPGTKCEIAG WGETKGTGND TVLNVALLN ISNQEENIKH RGHVRESEMC TEGLLAPVGA CEGDYGGPLA CFTHNCWVLK GIRIPNRVCT RSRWPAVFTR VSVFVDWIHK VMRLG
Source	Mammalian cell
Target Names	MST1L
Protein Names	Recommended name: Putative macrophage-stimulating protein MSTP9 Alternative name(s): Brain rescue factor 1 Short name= BRF-1 Hepatocyte growth factor-like protein homolog
Expression Region	21-715
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.