



Recombinant Human Uncharacterized protein KIAA1530 (KIAA1530)

Product Code	CSB-BP650553HU
Abbreviation	KIAA1530
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2YD98
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDQKLSKLV E LTTSGEPRL NPEKMKELKK ICKSSEEQLS RAYRLIIAQL TQEHAEIRLS AFQIVEELFV RSHQFRMLVV SNFQEFLELT LGTDPAQPLP PPREAAQRLR QATTRAVEGW NEKFG EAYKK LALGYHFLRH NKKVDFQDTN ARSLAERKRE EEKQKHLDKI YQERASQAER EMQEMSGEIE SCLTEVESCF RLLVPFDFDP NPETESLGMA SGMSDALRSS CAGQVGPCRS GTPDPRDGEQ PCCSRDLPAS AGHPRAGGGA QPSQTATGDP SDEDESDLE EFVRSHGLGS HKYTLDELVC SEGLKVQENE DNLALIHAAR DTLKLIRNKF LPAVCSWIQR FTRVGTGGC LKRAIDLKAE LELVLRKYKE LDIEPEGGER RRTEALGDAE EDEDEDVE VPEKEGYEPH IPDHLRPEYG LEAAPEKDTV VRCLRTRTRM DEEVSDPTSA AAQLRQLRDH LPPPSSASPS RALPEPQEAQ KLAERARAP VVPYGVDLHY WGQELPTAGK IVKSDSQHRF WKPSEVEEEV VNADISEMLR SRHITFAGKF EPVQHWCRAP RPDGRLCERQ DRLKCPFHGK IPRDDEGRP LPEDRAREQ RRQLQKQERP EWQDPELMRD VEAATGQDLG SSRYSGKGRG KKRRYPSLTN LKAQADTARA RIGRKVFAKA AVRRVVAAMN RMDQKKHEKF SNQFNIALN
Source	Baculovirus
Target Names	UVSSA
Protein Names	Recommended name: Uncharacterized protein KIAA1530
Expression Region	1-709
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.