



# Recombinant *Oryza sativa* subsp. japonica

## Probable protein phosphatase 2C 73

### (Os11g0109000, LOC\_Os11g01790)

<b>Product Code</b>	CSB-BP650870OFG
<b>Abbreviation</b>	Os11g0109000, LOC_Os11g01790
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q2RBJ6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGICCSKGKE ELEEEGFPWK HDAFFHDQLW SAGVSMHTKQ GWKGANQDAM TTCQDFAGHK GQIFCGVFDG HGPLGREVAR HVRDVLVVKL SSSLALKTEQ DPSSNTDKET LEKSDCTSLS DTSNEKQLLS TWKNIFVKTF EDVDEDLRQH SGIDCICSGT TAVTVVRQGD HLIANLGDS RAVLCTRDSK DRPISVQLTT DLKPNLPSEA ERILNSKGRV FAMDDEPDVP RMWLDPQDAP GLAMARAFGD FCLKSHGLIC TPEVYYRKLS AKDDFLVLAT DGIWDVLSNK EVIKIVSSAT DHSKAAKQLV ERAVRTWRRK FPTSMVDDCA VVCLFLKPSP SSESTPGDA KPPQAVSFTG SFRKVLGGGG GEAEEGTNVW RALEGVARVN SVVRLPRMGA VLSWRRRSTS LEEDDEARID
<b>Source</b>	Baculovirus
<b>Target Names</b>	Os11g0109000
<b>Protein Names</b>	Recommended name: Probable protein phosphatase 2C 73 Short name= OsPP2C73 EC= 3.1.3.16
<b>Expression Region</b>	1-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.