



Recombinant Human Translin-associated factor X-interacting protein 1 (TSNAXIP1)

Product Code	CSB-EP651118HU-B
Abbreviation	TSNAXIP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2TAA8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGGHLSPWPT YTSGQTILQN RKPCSDDYRK RVGSCQQHPF RTAKPQYLEE LENYLKELL LLDLGT DSTQ ELRLQPYREI FEFIEDFKT YKPLLSSIKN AYEGMLAHQR EKIRALEPLK AKLVTVNEDC NERILAMRAE EKYEISLLKK EKMNLLKLID KKNEEKISLQ SEVTKLRKNL AEEYLHYLSE RDACKILIAD LNELRYQRED MSLAQSPGIW GEDPVKLT LA LKMTRQDLTR TQMELNNMKA NFGDVVPRRD FEMQEKT N KD LQEQLDTLRA SYEEVRKEHE ILMQLHMSTL KERDQFFSEL QEIQRTSTPR PDWTKCKDVV AGGPERWQML AEGKNSDQLV DVLLEEIGSG LLREKDFFPG LGYGEAIPAF LRFDGLVENK KPSKKDVVNL LKDAWKERLA EEQKETFPDF FNFLEHRFG PSDAMAWAYT IFENIKIFHS NEVMSQFYAV LMGKRSENVY VTQKETVAQL LKEMTNADSQ NEGLLTMEQF NTVLKSTFPL KTEEIQIQLM EAGGWHPSSS NADLLNYRSL FMEDEEGQSE PFVQKLWEQY MDEKDEYLQQ LKQELGIELH EEVTLPLKLRG GLMTIDPSLD KQTVNTYMSQ AFQLPESEMP EEGDEKEEAV VEILQTALER LQVIDIRRVG PREPEPAS
Source	E.coli
Target Names	TSNAXIP1
Protein Names	Recommended name: Translin-associated factor X-interacting protein 1 Short name= Trax-interacting protein 1
Expression Region	1-658
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.