



Recombinant Solanum tuberosum Acetyl-coenzyme A carboxylase carboxyl transferase subunit beta, chloroplastic (accD)

Product Code	CSB-BP651219FIG
Abbreviation	accD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2VEG8
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	MERWWFNSML FKKEFERRCG LNKSMGSLGP IENTSEDPNL KMKNIHSCSN VDYLFGVKDI WNFISDDTFL VSDRNGDSYS IYFDIENQIF EVDNDHSFSL ELESSFSSYR NSSYLNNGFR GEDPYNSYM SYMYDTQYSW NNHINSCIDN YLQSQICIDT SIISGSESYG DSIYIRAICS GESLNSSENE GSSRRTRTKG SDLTIRESSN DLEVTQKYKH LWVQCENCYG LNYKKFLKSK MNICEQCGYH LKMSSSDRIE LLIDPGTWDP MDEDMVSLDP IEFHSEEEPY KDRIDSYQRK TGLTEAVQTG IGQLNGIPVA IGVMDFQFMG GSMGSVVGEK ITRLIEHAAN QNLPLIIVCA SGGARMQEGS LSLMQMAKIS SALYDYQLNK KLFYVSILTS PTTGGVTASF GMLGDIIIIE PNAYIAFAGK RVIEQTLNKT VPEGSQAAEY LFQKGLFDLI VPRNLLKSVL SELFKLHAFF PLNQKSSKIK
Source	Baculovirus
Target Names	accD
Protein Names	Recommended name: Acetyl-coenzyme A carboxylase carboxyl transferase subunit beta, chloroplastic Short name= ACCase subunit beta Short name= Acetyl-CoA carboxylase carboxyltransferase subunit beta EC= 6.4.1.2
Expression Region	1-490
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.