



# Recombinant *Oryza sativa* subsp. *japonica* CBL-interacting protein kinase 4 (CIPK4)

<b>Product Code</b>	CSB-MP651638OFG
<b>Abbreviation</b>	CIPK4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q2QMI0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAMAAMDARK KKS GGGGGEP LLGKYELGRM LGRGTFKVVY LARAVAGGEA VAVKVIDKAE VMGTAGMAPR VLREVAAMRR LRHPHVLRRLH EVLATRARIY LVMELATGGD LLSRLAALPR RRLPESAARR VVQLVDALS YCHARGVAHR DVKPQNVLLD GDGNLKVSDG GLAALPDTLR DDGRLHTACG TPAYAAPEVL RRRAYDGAKA DAWSCGVILF VLLAGHLPFD DSNIADMCRK AHRREYELPR WVSQPARRLV SRLDPNPDT RVAVESLAAH HPWFKRSLSV DSQLDGLLNG EPERAVAFQA APPPPLNAFD IISMSPGLDL SGLFGEHDKS LREKRFTTTA SPEKTLEQLG LAGGKLGYYV VVGKKGVECL PLAGGRLSSG IAAMSVEMSE VAPPLLLVEL RLEVAAGDVD GGDGEVKGFG WEQLRMELGD VVRAWHSCED LCEI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CIPK4
<b>Protein Names</b>	Recommended name: CBL-interacting protein kinase 4 EC= 2.7.11.1 Alternative name(s): OsCIPK04
<b>Expression Region</b>	1-454
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.