



Recombinant *Oryza sativa* subsp. japonica Putative germin-like protein 12-4 (Os12g0155000, LOC_Os12g05880)

Product Code	CSB-EP651645OFG
Abbreviation	Os12g0155000, LOC_Os12g05880
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2QXJ0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	SDPSPLQD FCVADKHSPV RVNGLPCKDA KDVSVDFFL AANLDKPMDT TKSKAGSNVT LINVMKLAGL NTLISMARI DYAPKGGQNP HTHPRATEIL TVLEGSLYVG FVTSNQANRE NKLFTKTLNK GDVVFVFPQGL IHFQFNPSYD KPAVAIAALS SQNPGAITIA NAVFGSHPI SDDVLAKAFQ VDKKAMDWLQ AQFWENNHN
Source	E.coli
Target Names	Os12g0155000
Protein Names	Recommended name: Putative germin-like protein 12-4
Expression Region	23-229
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.