



Recombinant *Oryza sativa* subsp. japonica PHD finger protein ALFIN-LIKE 8 (Os11g0244800, LOC_Os11g14010)

Product Code	CSB-MP651652OFG
Abbreviation	Os11g0244800, LOC_Os11g14010
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2R837
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDGGGAHRTP EDVFRDFRAR RAGMIKALTT DVEKFYQQCD PEKENLCLYG LPNETWEVNL PAEEVPELP EPALGINFAR DGMDEKDWLS LVAVHSDTWL LAVAFYFGAR FGFDKESRKR LFSMINNLT IYEVVTGTAK QKSKEKTPKT SGKSNKSGTK PSRQPEPNSR GPKMPPPKDE DDSGGEEEE EEDHENTLCG ACGDNYGQDE FWICCDACET WFHGKCVKIT PAKAEHIKHY KCPNCSSSSK RARA
Source	Mammalian cell
Target Names	Os11g0244800
Protein Names	Recommended name: PHD finger protein ALFIN-LIKE 8
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.