



Recombinant Human Brorin (VWC2)

Product Code	CSB-BP651907HU
Abbreviation	VWC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2TAL6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SPS IPLEKLAQAP EQPGQEKREH ASRDGPGRVN ELGRPARDEG GSGRDWKS KS GRGLAGREP W SKLKQAWVSQ GGGAKAGDLQ VRPRGDT PQA EALAAAQDA IGPELAPTPE PPEEYVYPDY RGKGC VDESG FVYAIG EKFA PGPSAC PCLC TEEGPLCAQP ECPRLHPRCI HVDTSQCCPQ CKERKN YCEF RGKTYQTLEE FVSPCERCR CEANGEV LCT VSACPQTECV DPVYEPDQCC PICKNGP NCF AETAVIPAGR EVKTDECTIC HCTYE EGTWR IERQAMCTR H ECRQM
Source	Baculovirus
Target Names	VWC2
Protein Names	Recommended name: Brorin Alternative name(s): Brain-specific chordin-like protein von Willebrand factor C domain-containing protein 2
Expression Region	28-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.