



Recombinant *Aspergillus oryzae* Autophagy-related protein 3 (atg3)

Product Code	CSB-BP651941DPA
Abbreviation	atg3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2U7R4
Product Type	Recombinant Protein
Immunogen Species	<i>Aspergillus oryzae</i> (strain ATCC 42149 / RIB 40) (Yellow koji mold)
Purity	>85% (SDS-PAGE)
Sequence	MNILHSTLST WRDRLAPVSR TSTFRNTGQI TPEEFVLAGD YLVYKFPSWS WADASNPAKR VSYLPPGKQF LVTRGVPCHR RLNDNFAGDA GHDELVRDM LSGGTGGVDD DGWLRTGGGQ DSADRQENRI KDVRTVDESG NMGEREEEEED EIPDMEDEDD DEEAIIRDPA SGTTPQTRTY NLYITYSNFY RTPRLYMSGY LSPSEPLPPH LMMEDVVG DY KDKVTLEDF PWYDGNVKMA SVHPCRHASV MKTLLDRADA ALKLRREKLK QAQSDPSKAP SVGESGLEGL VDDIKALSLS DQQQHGS DKS GGDEWEVLQH DEEEQVAIRV DQYLVVFLKF IASVTPGIEH DFTMGV
Source	Baculovirus
Target Names	atg3
Protein Names	Recommended name: Autophagy-related protein 3 Alternative name(s): Autophagy-related E2-like conjugation enzyme atg3
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.