



Recombinant Human Leucine-rich repeat-containing protein FAM211B (FAM211B)

Product Code	CSB-BP652020HU
Abbreviation	FAM211B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2VPJ9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGARLGRRAG PEAGSEAGAA AGCGPAPYER RVRWLREIQS TLRERRPERA RQLLRLLRQD LGLERTLLPD ILYRDVAFLN PVDPIHDL VNLARDLQCP KKDYELWKSS DKICRQLIYH LTPHSKQQQG SSLRQRKTQS CLKSSLQKTL LAGETVDLSG IPLSTQDVQH ITRYLSSHGA VLAVLDLSFT GLSDELLHLL LPSLWALPRL TQLLLNGNRL TRATARKLTD AIKDTTKFPA LAWVDLGNNV DVASLPQPLL VGLRRRLSQR TSLPTIYEGL DLEPEGSAAG ATTPASTWDS TAAGLGPEPQ ACCAR
Source	Baculovirus
Target Names	LRRC75B
Protein Names	Recommended name: Leucine-rich repeat-containing protein FAM211B
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.