



Recombinant Human Kelch-like protein 38 (KLHL38)

Product Code	CSB-MP652063HU
Abbreviation	KLHL38
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2WGJ6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDEESLDGLL FKDHDFSSDL LRQLNSLRQS RILTDVSICA GAREIPCHRN VLASSSPYFR AMFCSSFREK SEAKVQLKGI DPPTLDQIVS YVYTGEAHIA TDNVLPVMEA ASMLQFPKLF EACSSYLQSQ LAPSNCLGMI RLSEILSCET LKKKAREVAL TSFPEVAASA DLKELCALEL RDYLGDDGLC GEEEKVFEAL MVWIKHDLQA RKRYMQELFK QVRLQYIHPA FFHHFIANDA LLQSSPACQI ILETAKRQMF SLCGTTVPDC KLLLVHPPRN SYQDFLILLG GRKDSQQTTR DVLLYSKQTG QWQSLAKLPT RLYKASAITL HRSIYVLGGM AVSSGRSLVS HNVIYFSLKL NQWRLGEPML VARYSHRSTA HKNFIFSIGG IGEGQELMGS MERYDSICNV WESMASMPVG VLHPAVAVKD QRLYLFGGED IMQNPVRLIQ VYHISRNSWF KMETRMKNV CAPAVVLGER IVIVGGYTRR ILAYDPQSNK FVKCADMKDR RMHHGATVMG NKLYVTGGRR LTTDCNIEDS ASFDCYDPET DTWTSQGQLP HKLFDHACL T LQCIPRTSGL P
Source	Mammalian cell
Target Names	KLHL38
Protein Names	Recommended name: Kelch-like protein 38
Expression Region	1-581
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.