



# Recombinant Bovine Origin recognition complex subunit 4 (ORC4)

<b>Product Code</b>	CSB-BP652122BO
<b>Abbreviation</b>	ORC4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q2YDI2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSSRKS KSSN LIQADYISQV QRILRERFCH QSPHGNLFGV QVQYKHLIEL LKRTAIHGES NSILIIGPRG SGKTM LINHA LKELMEIEGV SENILQVHLN GLLQINDKIA LKEITRQLNL ENVVGD KVFVG SFAENLSFLL EALKKGDRTS SCPVIFILDE FDLFAHHKNQ TLLYNLLDIS QSAQTPVVII GLTCRLDILE LLEKRVKSRF SHRQIHLMNS FSFPQYLKIF KEQLSLPSVF PEEIFAEKWN ENVQSLSEDR SVKEVLQKHF NVSKNLRSLH MLLMLALNRV TPSHPFITAA DLMEANQLCS MDSKGNIVHG LSVLEICLII AMKHLNDIYE EEPFNFQMVY NEFQKFVQRK AHSVYNFEKP VVMKAFEHLQ QLELIRPMER TSVNAQREYQ LMKLLLDNTQ IMNALQKYPN CPTDVRQWAT SSSLWL
<b>Source</b>	Baculovirus
<b>Target Names</b>	ORC4
<b>Protein Names</b>	Recommended name: Origin recognition complex subunit 4
<b>Expression Region</b>	1-436
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.