



Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 66 (Os12g0405100, LOC_Os12g21700)

Product Code	CSB-BP652423OFG
Abbreviation	Os12g0405100, LOC_Os12g21700
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2QT65
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAAGAGAGGG GEGDSNGGG TSPGGVSAAA PAIGPHHLGV AAAEEMWQM TLGGGESMES TPYPERIGEP DCSYYMRTGL CRFGMTCKFN HPPNRKLAVA AARMNGEYPY RVGQPECQYY LKTGTCKFGA TCKFHHPREK AALANRVQLN VLGYPMRPNE KECAYYLRTG QCKFASTCKF HHPQPSNTMV AVRNSMYSYG QSATSPGQHT YPGAVTNWTL SRSASFIASP RWPGHSGYQAQ VIVPQGLVQV PGWNPYAAQM GSSSPDDQQR TPVTTQYYGS RQSETGGMGD HGMYQSYQGG SVPVGVYTVQ GENIFPERPD QPECQFYMKT GDCKFGAVCK FHHPKERLVP APNCALNSLG LPLRPGEPVC TFYSRYGICK FGPNCKFDHP MGTLMYGSAT SPRGDVSSMH YQLSPSPGHP GILLDGGSGR SHRVPQSDSQ QIPSGDGNAE REAS
Source	Baculovirus
Target Names	Os12g0405100
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 66 Short name= OsC3H66
Expression Region	1-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.