



Recombinant *Oryza sativa* subsp. japonica Probable protein phosphatase 2C 75 (Os11g0417400, LOC_Os11g22404)

Product Code	CSB-EP652433OFG
Abbreviation	Os11g0417400, LOC_Os11g22404
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2R637
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGTCLTTAEQ RAMEVPAASV KGGGRRSDE EAPGRIAGNG AGNVACLFTR QGKKG TNQDA MVAWENYNGR SDTVFCGVFD GHGPHGH LIA RKVRDILPSR LCDLIYEDCG DSPTSNSDVS TLEENLSPYA DAECRSPTLA GQKEHQEFFN AMKESFRKAF KNVDKELKLQ RNIDSICSGT TAVTLIKQGH DLIVGNLGD S RAVLGTRDQN DKLVAHQ LTV DLKPDHPREA RRIRRCNGRV FAHQDEPDVA RLWLPNCNSP GLAMARAFGD FCLKDFGLIS VPDVTYRQIT EKDEFIVLAT DGVWDVLSNQ EVVDVVASCS GRFAAARSVV DLANETWRFK YPTSKTDDCA VVCLFLNKYE VTGGLSGQPG YSPRMPALSG ITRPNSKRVT PDDVDDGSDS NVSGDERSLD GFTRLN TLLA LPKFGDTSPT KK
Source	E.coli
Target Names	Os11g0417400
Protein Names	Recommended name: Probable protein phosphatase 2C 75 Short name= OsPP2C75 EC= 3.1.3.16
Expression Region	1-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.