



Recombinant Human Prickle-like protein 4 (PRICKLE4)

Product Code	CSB-MP652689HU
Abbreviation	PRICKLE4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2TBC4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSPQGPAVLS LGSLCLDTNQ APNWTGLQTL LQQLPPQDID ERYCLALGEE ERAELQLFCA RRKQEALGQG VARLVLPKLE GHTCEKCREL LKPGEYGVFA ARAGEQRCWH QPCFACQACG QALINLIYFY HDGQLYCGRH HAELLRPRCP ACDQLIFSWR CTEAEGQRWH ENHFCCQDCA GPLGGGRYAL PGGSPCCPSC FENRYSDAGS SWAGALEGQA FLGETGLDRT EGRDQTSVNS ATLSRTLLAA AGGSSLQTQR GLPGSSPQQE NRPGDKAEAP KGQEQCRLET IRDPKDTPFS TCSSSSDSEP EGFFLGERLP QSWKTPGSLQ AEDSNASKTH CTMC
Source	Mammalian cell
Target Names	PRICKLE4
Protein Names	Recommended name: Prickle-like protein 4 Alternative name(s): Overexpressed breast tumor protein
Expression Region	1-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.