



# Recombinant Arabidopsis thaliana Uncharacterized zinc finger protein At4g06634 (At4g06634)

<b>Product Code</b>	CSB-BP652768DOA
<b>Abbreviation</b>	At4g06634
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q2V3L3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDHQNYQYQN PFERRPILKS KAPAVKWIKE WVPQDIVATG GKCHLHKWVT EDTFSRLKEK EKEPDVPEPE PEPTTEILFL CSYDGCCKTF FDVSALRKHS HIHGERQYVC DQEGCGKKFL DSSKLRHYL IHTGERNYIC TYEGCGKAFS LDFNLRSHMK THSQENYHIC PYSGCVKRYA HEYKLNHVA AYHEKNGGGE TPKYTPPAEK VLRTVKTPAT VCGPSSDRPY ACPYEGCEKA YIHEYKLLKH LKREHPGHLQ EENADTPTLN KHNGNDRNEI DDGSDQDVYR KHASNGKGQT HKQQSRAKPN MRTPPAKVGK KGSTSSPAKA RIAKKPWQAK ETFEEVEREE EEDSEETEED RDNVEDGWRF GENNEDDDDD EETEYED
<b>Source</b>	Baculovirus
<b>Target Names</b>	YY1
<b>Protein Names</b>	Recommended name: Uncharacterized zinc finger protein At4g06634
<b>Expression Region</b>	1-387
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.