



# Recombinant Danio rerio Threonine synthase-like 2 (thnsl2)

<b>Product Code</b>	CSB-BP652902DIL
<b>Abbreviation</b>	thnsl2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q2YDP8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRYCSTRAGV QGRSFRDVLV SGYAADGGMF MPENLPSLSA ETLSWRRLS YRQLLCEVCE LFIPEQEIPR EQLEGLVSSA LSSFVSPDAV RLVQLKDSLS ILELFHGETL AFKDLAMSCCT AHFLQYFLRR DRQRATILVG TSGDTGSSAI RSVLGLREVD IVVVFPGRRI TKIQELQMTT SVAENVHVFA ADGTSDDDIDV PLRKLAFADAD LVQRHRLMSL NSVNWSRIMV QTAHFLFAYL QLTPSLPEGD TLPVLEVLVP TGGAGNITAG IIVKRMGVPL RLVAMVNAND IVHRTVQSGD FSMSSSVTQT LAPAIDIQDP YNMERVFWLL SGGDGLMVKS LMEEFQKTHK LSLPASLHQQ LSEVLSSGSV SDDGILEAMR RCWQDNHYLI CPHTAVAVWR HYQSPVRPGE SRCCIATASP HKFQQAVQRA GLTLELPESL QVLDQLETRC ADLKCSDDWE ETLRRHIESI SSRNS
<b>Source</b>	Baculovirus
<b>Target Names</b>	thnsl2
<b>Protein Names</b>	Recommended name: Threonine synthase-like 2 Short name= TSH2 EC= 4.2.3.-
<b>Expression Region</b>	1-476
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.