



# Recombinant *Synechococcus elongatus* FO synthase subunit 2 (cofH)

<b>Product Code</b>	CSB-EP653236FPY
<b>Abbreviation</b>	cofH
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q31QF1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Synechococcus elongatus</i> (strain PCC 7942) ( <i>Anacystis nidulans</i> R2)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MIDSATVAAI LASVLDGKPL EPEAATVLLK ARDRSLRQQI QAAANQLRSR QVGDRVSYVI NRNLNFTNIC EQHCNFCFR RDADQDGAFW LDASILLEKG AAAVAAGATE FCLQGGLNPA AKRNGRSLDF YVELTASLKQ AFPQIHLHAF SPQEIQFIAR EDGLSFREVL MALRSAGVGS LPGTAAEVLD DSVRRILCPE KLDSATWKTI IQTAHQVGLP TTSTLLSGHL ETPSQQAQHL EQLRQLQAAA IAGETPARIT EFILLPFVGE LAPAPLRKRV KRDQPDLSA LLVMAVARLY LGDWIANHQP SWVKLGLAGA TQALDWGCND LGGTLMEEHI TSMAGAQQGT AQTVEQLEAA IAAAGRQPYQ RDTLYRPVAV EAVHAG
<b>Source</b>	E.coli
<b>Target Names</b>	cofH
<b>Protein Names</b>	Recommended name: FO synthase subunit 2 EC= 2.5.1.77 Alternative name(s): 7,8-didemethyl-8-hydroxy-5-deazariboflavin synthase subunit 2
<b>Expression Region</b>	1-386
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.