



Recombinant *Synechococcus elongatus* FO synthase subunit 2 (cofH)

Product Code	CSB-YP653236FPY
Abbreviation	cofH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31QF1
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MIDSATVAAI LASVLDGKPL EPEAATVLLK ARDRSLRQQI QAAANQLRSR QVGDRVSYVI NRNLNFTNIC EQHCNFCFR RDADQDGAFW LDASILLEKG AAVAAGATE FCLQGGLNPA AKRNGRSLDF YVELTASLKQ AFPQIHLHAF SPQEIQFIAR EDGLSFREVL MALRSAGVGS LPGTAAEVLD DSVRRILCPE KLDSATWKTIIQTAHQVGLP TTSTLLSGHL ETPSQQAQHL EQLRQLQQA IAGETPARIT EFILLPFVGE LAPAPLRKRV KRDQPDLSA LLVMAVARLY LGDWIANHQP SWVKLGLAGA TQALDWGCND LGGTLMEEHI TSMAGAQQGT AQTVEQLEAA IAAAGRQPYQ RDTLYRPVAV EAVHAG
Source	Yeast
Target Names	cofH
Protein Names	Recommended name: FO synthase subunit 2 EC= 2.5.1.77 Alternative name(s): 7,8-didemethyl-8-hydroxy-5-deazariboflavin synthase subunit 2
Expression Region	1-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.