



Recombinant *Synechococcus elongatus* Thiazole synthase (thiG)

Product Code	CSB-BP653238FPY
Abbreviation	thiG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31QR1
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	≥85% (SDS-PAGE)
Sequence	MENSCADAEF MVATPDDTTAD VLTIAGRSFR SRLMTGTGKY RSFEQMRASI AASGCEIVTV AVRRVQTNAP GHEGLAEALD WQKIWMLPNT AGCATAEEAI RVARLGREMA KLLGQEDNNF VKLEVIPDSR YLLPDPIGTL QAAEQLVKEG FAVLPIYINAD PLLAKRLEEY GCATVMPLGS PIGSGQGIRN EANIRIIVEN AKVPVVVDAG IGTASEAAQA MELGADALLI NTAIAQAEDP ARMAQAMAMA TIAGRLAFQA GRIPTRSAAI ASSPQTGLVG QSPATV
Source	Baculovirus
Target Names	thiG
Protein Names	Recommended name: Thiazole synthase EC= 2.8.1.10
Expression Region	1-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.