



# Recombinant Saxifraga cernua Maturase K (matK)

<b>Product Code</b>	CSB-MP653520SCAL
<b>Abbreviation</b>	matK
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q33078
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Saxifraga cernua (Nodding saxifrage) (Lobaria cernua)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGEFRRYFEL DKFLQQDFLY PLIFHEYIYA LAHDHVLNRF ILVDNFGYDN KSSSLMVKRL ITRMSQQNPL LISANDSNXK KWLGHKNKNLY SQMISEGFTV IGEIPFSLRL VSSVERKEEI GKSHNLSIH SLFPFLEDNF VHLNYVSDIL IPYPIHLEIV VQTLRSWLKD ASSHLRLRFL LYEQNWNSL IIPTPKKPSS IVSKRNQRLF LFLYNSYVCE YESIFFFLCN QSSHLPSTSS ETLFERTFFH GKIKALVEVG SNDFFPVRWL FKDPFMHYVR YQGKSILASK STPLLMDKWK YYLVNFWQCH LYVWSKPVRI HMNQLSNHSF GFLGYLSSVR LNPSVVRSQM LE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	matK
<b>Protein Names</b>	Recommended name: Maturase K Alternative name(s): Intron maturase
<b>Expression Region</b>	1-352
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.