



Recombinant Arabidopsis thaliana PAP-specific phosphatase HAL2-like (AHL)

Product Code	CSB-MP653633DOA
Abbreviation	AHL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q38945
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAVDSLETEI DTAVRVVHLA SSLCVKVQEK LHLPNGGHVK SKDDDSPVTV ADFGVQAIVS WVLAEVFGDQ NLSIVAEEDT ETLSEADSLG LLGAVSNAVN EALSEAQNYG LPKPKPLGS SEILKAISRC NSVGGPKGRH WVLDPVDGTL GFVRGDQYAV ALALIENGKV LLGVLGCPNY PVKKECLSNG CNQAMKTKAV AGSVSKGCVM YAKRGSGQAW MQPLIVGGIP ESATLLKVSS VDDPVLATVC EPVERANSNH LFTAGLANSM GVRKQPMRVY SMVKYAAIAR GDAEVFMKFA QSSYKEKIWD HAAGVVIVEE AGGVVTDAGG RNLDFSKGVY LEGLDRGIIA CSGQVLHEKI IGAVYASWES SSL
Source	Mammalian cell
Target Names	AHL
Protein Names	Recommended name: PAP-specific phosphatase HAL2-like Alternative name(s): 3'(2'),5'-bisphosphate nucleotidase EC= 3.1.3.7 3'(2'),5'-bisphosphonucleoside 3'(2')-phosphohydrolase DPNPase Halotolerance protein
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.