



Recombinant *Synechococcus elongatus* UPF0176 protein Synpcc7942_1892 (Synpcc7942_1892)

Product Code	CSB-EP653925FPY
Abbreviation	Synpcc7942_1892
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31LZ7
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MSLVLINFYR FVALGDCDRW RQWLQDLCTA LGLRGTILLA PEGINAGLAG NTEAIAQFLS ELQQHPPFAN LSFKSATVTD WPFARLKVKV KPEIVSLGCP ELNPAERTGT LVAPQDWNQL LQDPEVVLID VRNRFEIALG SFPRAIDPQT DRFRDFPRFV QEQLLPQPPA KVAMFCTGGI RCEKASAYLL EQGIETVYQL EGGILNYLEA IAPEENHWQG DCFVFDERIA VDRQLQTPQH QLCPCAGQPV VATTCSHCQD SVQASSSPK
Source	E.coli
Target Names	Synpcc7942_1892
Protein Names	Recommended name: UPF0176 protein Synpcc7942_1892
Expression Region	1-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.