



Recombinant *Synechococcus elongatus* Ribulose biphosphate carboxylase large chain (cbbL)

Product Code	CSB-MP653931FPY
Abbreviation	cbbL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31NB3
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MPKTQSAAGY KAGVKDYKLT YYTPDYTPKD TDLLAAFRFS PQPGVPADEA GAAIAAESST GTWTTVWTDL LTDMDRYKKGK CYHIEPVQGE ENSYFAFIAY PLDLFEEGSV TNILTSIVGN VFGFKAIRSL RLEDIRFPVA LVKTFQGPPH GIQVERDLLN KYGRPMLGCT IKPKLGLSAK NYGRAVYECL RGGLDFTKDD ENINSQPFQR WRDRFLFVAD AIHKSQAETG EIKGHYLNVT APTCEEMMKR AEFAKELGMP IIMHDFLTAG FTANTTLAKW CRDNGVLLHI HRAMHAVIDR QRNHGIIHFRV LAKCLRLSGG DHLHSGTVVG KLEGDKASTL GFVDLMREDH IEADRSRGVF FTQDWASMPG VLPVASGGIH VWHMPALVEI FGDDSVLQFG GGTLGHPWGN APGATANRVA LEACVQARNE GRDLYREGGD ILREAGKWSP ELAAALDLWK EIKFEFETMD KL
Source	Mammalian cell
Target Names	cbbL
Protein Names	Recommended name: Ribulose biphosphate carboxylase large chain Short name= RuBisCO large subunit EC= 4.1.1.39
Expression Region	1-472
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.