



Recombinant *Synechococcus elongatus* Phosphoribosylglycinamide formyltransferase 2 (purT)

Product Code	CSB-MP653941FPY
Abbreviation	purT
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31QP9
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MTFADRLPRR LMLLGSSELG KEFAIAAQRL GNTVIAVDY AHAPAMQVAD ASAVISMLDG EALEAVVQEF QPDLIPEIE AIRTEKLQEF EDRGLTVIPT ARATHFTMNR DRIRDLAAQQ LGLRTARYAY ASCFEELQTV AAAIGYPNVI KPVMSSSGKG QSIIQQPDQL QAAWDYAIAG SRGDSQKVIL EEFIPFELEI TLLTIRQWQG PTLFCPPIGH RQERGDYQES WQPAPLRPDL LAQAQIAAQ VTEALGGAGL FGVEFFVTPD EVIFSELSPR PHDTGMVTLI SQNLNEFELH LRAVLGLPIP AIELLGPSAS RVILAEDSGD RPSYAGVAAA LQEPWVDLRL FGKPDMPQR RMGVALARDE SVEAARAKAD RAASQVTIQP S
Source	Mammalian cell
Target Names	purT
Protein Names	Recommended name: Phosphoribosylglycinamide formyltransferase 2 Short name= GART 2 EC= 2.1.2.- Alternative name(s): 5'-phosphoribosylglycinamide transformylase 2 Formate-dependent GAR transformylase GAR transformylase 2
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.