



Recombinant Bovine Nucleotide exchange factor SIL1 (SIL1)

Product Code	CSB-EP654161BO
Abbreviation	SIL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q32KV6
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	HKNSEFALT NPEKSSTKET ERKETEEEELE PKILEVFHPT HEWQALRPGQ AVPAGSHVRL NLQTGAREVK LHDEDKFQTN LKGLKKGKRL DINTNTYTSQ DLKNALAKFK EGAEMENSKE DQERQAKVKR LFRPIEELKK DFEELNVVIE TDMQIMVRLI NKFNSSSSSL EEKIAALFDL EYYVHQMDNA QDLLSFGGLQ VVINGLNSTE PLVKEYAAFV LGAAFSSNPK VQVEAIEGGA LQKLLVILAT EQPLTTKKKA LFALCSLLRH FPYAQQQFLK LGGLQVLRSL VQEKGMVLA VRVVTLLYDL VTEKMFAEED AELTWDMSPK RLQQYRQVHL LPSLQEQGWC EITAHLLALP EHDAREKVLQ TLGALLATCR DRFHQDTQLH RTLGLTQAEY QALAALELQE GEDEGYFREL LDSIDSLLTE LR
Source	E.coli
Target Names	SIL1
Protein Names	Recommended name: Nucleotide exchange factor SIL1
Expression Region	32-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a resident endoplasmic reticulum (ER), N-linked glycoprotein with an N-terminal ER targeting sequence, 2 putative N-glycosylation sites, and a C-terminal ER retention signal. This protein functions as a nucleotide exchange factor for another unfolded protein response protein. Mutations in this gene have been associated with Marinesco-Sjogren syndrome. Alternate transcriptional splice variants have been characterized.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.