



Recombinant Arabidopsis thaliana Rac-like GTP-binding protein ARAC4 (ARAC4)

Product Code	CSB-EP654338DOA-B
Abbreviation	ARAC4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q38919
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MASRFIKCVT VGDGAVGKTC MLISYTSNTF PTDYVPTVFD NFSANVVVDG NTVNLGLWDT AGQEDYNRLR PLSYRGADV ILAFSLISKA SYENIAKKWI PELRHYAPGV PIILVGTKLD LRDDKQFFID HPGAVPITTN QGEELKKLIG SAVYIECSSK TQQNVKAVFD AAIKVVLQPP KQKKKKKNKN RC
Source	E.coli
Target Names	ARAC4
Protein Names	Recommended name: Rac-like GTP-binding protein ARAC4 Alternative name(s): GTPase protein ROP2
Expression Region	1-192
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.