



Recombinant Arabidopsis thaliana Shaggy-related protein kinase iota (ASK9)

Product Code	CSB-BP654391DOA
Abbreviation	ASK9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q39012
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MASLPLGPQP HALAPPLQLH DGDALKRRPE LDSDKEMSAA VIEGNDAVTG HIISTTIGGK NGEPKQTISY MAERVVGTGS FGIVFQAKCL ETGESVAIKK VLQDRRYKNR ELQLMRPMDH PNVISLKHCF FSTTSRDELFLNLVMEYVPE TLYRVLRYHT SSNQRMPIFY VKLYTYQIFR GLAYIHTVPG VCHRDKVPQN LLVDPLTHQV KLCDFGSAKV LVKGEPNISY ICSRYRRAPE LIFGATEYTA SIDIWSAGCV LAELLLGQPL FPGENSVDQL VEIIKVLGTP TREEIRCMNP NYTDFRFPQI KAHPWHKVFH KRMPPEIDL ASRLLQYSPS LRCTALEACA HPFFNELREP NARLPNGRPL PPLFNFKQEL GGASMELINR LIPEHVRQRM STGLQNS
Source	Baculovirus
Target Names	ASK9
Protein Names	Recommended name: Shaggy-related protein kinase iota EC= 2.7.11.1 Alternative name(s): ASK-iota
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.