



# Recombinant *Synechococcus elongatus* Glutamyl-tRNA (Gln) amidotransferase subunit A (gatA)

<b>Product Code</b>	CSB-MP654627FPY
<b>Abbreviation</b>	gatA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q31LC2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Synechococcus elongatus</i> (strain PCC 7942) ( <i>Anacystis nidulans</i> R2)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASIRELHSQ LVRKERSAVE IAEAALQRIE TLEPQLKSFL HVTADQAIQ AKAVDQRIAA GEEISLLAGI PIGIKDNLCT RGIPTTCASK ILQGFVPPYE STVTQRLAEA GAVAVGKTNL DEFAMGSSTE NSAYQTTANP WDLTRVPGGS SGGSAAAVAA DECVVALGSD TGG SIRQPAA FCGVVGLKPT YGLVSRVGLV AYASSLDQIG PFSRSVEDTA ILLGAIAGHD PKDATSLKVE VPDYTQFLKP SLAGIKVGV I TETVTDSPAG QAMQAALV L QGLGAEIREI SCPRFAYGLP AYYIIAPSEA SANLARYDGV KYGYRVEEAE TLIDMYCRTR AEGFGSEVKR RIMIGTYALS AGYYDAYYLK AQKVRTLIQ DFEEAFAEVD VLVSPPTPTT AFKAGEKTAD PLSMYLSDLM TIPVNLGLP GLSLPCGFDE AGLPYGLQII GNVLRDQVL HTAYAYEQAT EWHLRQPAL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	gatA
<b>Protein Names</b>	Recommended name: Glutamyl-tRNA(Gln) amidotransferase subunit A Short name= Glu-ADT subunit A EC= 6.3.5.-
<b>Expression Region</b>	1-479
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.