



Recombinant *Synechococcus elongatus* Photosystem II reaction center Psb28 protein (psb28)

Product Code	CSB-MP654633FPY
Abbreviation	psb28
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31ML0
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MSAAIQFTRG IDEPVVPDVR LTRSRDGQNG TATFYFDEPQ ALVGEVRQDI TGMYLDDDEG ELATREVKAK FINGQPAGLE AVYLMRSPEE WDRFMRFMQR YAEANGLGFT EA
Source	Mammalian cell
Target Names	psb28
Protein Names	Recommended name: Photosystem II reaction center Psb28 protein Alternative name(s): Photosystem II 13 kDa protein Photosystem II reaction center W protein
Expression Region	1-112
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.