



Recombinant *Synechococcus elongatus* Tryptophan synthase alpha chain (trpA)

Product Code	CSB-MP654646FPY
Abbreviation	trpA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31R76
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	≥85% (SDS-PAGE)
Sequence	MTAISDCFAA LRSQGRCALI PFLTAGDPDL ETTRQALLAL DREGADLIEL GVPYSDPLAD GPVIQAAATR ALQAGTRLDD VLALLKDVRS QIKAPIVLFT YCNPILNRGF EAFLDQIAAA GANGLVVPDL PLEESQRLSE VAAERGIDLI LLIAPTSSAD RIAAISKQAR GFIYLVSVTG VTGMRQGMQS RVADLLQEIR QGTDKPIGVG FGISGAEQAR QVRDWGADGV IVGSAFVNRL QEQGVGVAT LCRELRQAID RQPVLS
Source	Mammalian cell
Target Names	trpA
Protein Names	Recommended name: Tryptophan synthase alpha chain EC= 4.2.1.20
Expression Region	1-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.