



Recombinant Shigella dysenteriae serotype 1 Lipoyl synthase (lipA)

Product Code	CSB-MP654846SGF
Abbreviation	lipA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q32IV2
Product Type	Recombinant Protein
Immunogen Species	Shigella dysenteriae serotype 1 (strain Sd197)
Purity	>85% (SDS-PAGE)
Sequence	MSKPIVMERG VKYRDADKMA LIPVKNVATE REALLRKPEW MKIKLPADST RIQGIKAAMR KNGLHSVCEE ASCPNLAECF NHGTATFMIL GAICTRRCPF CDVAHGRPVA PDANEPVKLA QTIADMALRY VVITSVDRDD LRDGGAQHFA DCITAIREKS PQIKIETLVP DFRGRMDRAL DILTATPPDV FNHNLENVPR IYRQVRPGAD YNWSLKLLER FKEAHPEIPT KSGLMVGLGE TNEEIIIVMR DLRRHGVTML TLGQYLQPSR HHLVPQRYVS PDEFDEMKAELAMGFTTHAA CGPFVRRSSYH ADLQAKGMEV K
Source	Mammalian cell
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.