



# Recombinant Human Zinc finger protein 3 (ZNF3)

<b>Product Code</b>	CSB-BP026651HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P17036
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	METQADLVSQ EPQALLDSAL PSKVPAFSDK DSLGDEMLAA ALLKAKSQEL VTFEDVAVYF IRKEWKRLEP AQRDLYRDVM LENYGNVFSL DRETRTENDQ EISEDTRSHG VLLGRFQKDI SQGLKFKEAY EREVSLKRPL GNSPGERLNR KMPDFGQVTV EEKLTTPRGER SEKYNDFGNS FTVNSNLISH QRLPVGDRPH KCDECSKSFN RTSDLIQHQR IHTGEKPYEC NECGKAFSQS SHLIQHQR IH TGEKPYECSD CGKTFSCSSA LILHRRIHG EKPYECNECG KTFWSWSTLT HHQRIHTGEK PYACNECGKA FSRSSLIHH QRIHTGEKPY ECNECGKA FS QSSHLYQHQR IHTGEKPYEC MECGGKFTYS SGLIQHQR IH TGENPYECSE CGKAFRYSSA LVRHQRIHTG EKPLNGIGMS KSSLRVTTTEL NIREST
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZNF3
<b>Protein Names</b>	Recommended name: Zinc finger protein 3 Alternative name(s): Zinc finger protein HF.12 Zinc finger protein HZF3.1 Zinc finger protein KOX25
<b>Expression Region</b>	1-446
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.