



Recombinant Human Zinc finger protein 3 (ZNF3)

Product Code	CSB-MP026651HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P17036
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	METQADLVSQ EPQALLDSAL PSKVPAFSDK DSLGDEMLAA ALLKAKSQEL VTFEDVAVYF IRKEWKRLEP AQRDLYRDVM LENYGNVFSL DRETRTENDQ EISEDTRSHG VLLGRFQKDI SQGLKFKEAY EREVSLKRPL GNSPGERLNR KMPDFGQVTV EEKLTTPRGER SEKYNDFGNS FTVNSNLISH QRLPVGDRPH KCDECSKSFN RTSDLIQHQR IHTGEKPYEC NECGKAFSQS SHLIQHQR IH TGEKPYECSD CGKTFSCSSA LILHRRIHG EKPYECNECG KTFWSWSTLT HHQRIHTGEK PYACNECGKA FSRSSLIHH QRIHTGEKPY ECNECGKA FS QSSHLYQHQR IHTGEKPYEC MECGGKFTYS SGLIQHQR IH TGENPYECSE CGKAFRYSSA LVRHQRIHTG EKPLNGIGMS KSSLRVTTTEL NIREST
Source	Mammalian cell
Target Names	ZNF3
Protein Names	Recommended name: Zinc finger protein 3 Alternative name(s): Zinc finger protein HF.12 Zinc finger protein HZF3.1 Zinc finger protein KOX25
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.