



Recombinant Arabidopsis thaliana Rac-like GTP-binding protein ARAC2 (ARAC2)

Product Code	CSB-YP655038DOA
Abbreviation	ARAC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q38903
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSTARFIKCV TVGDGAVGKT CMLISYTSNT FPTDYVPTVF DNFSANVVVD GSTVNLGLWD TAGQEDYNRL RPLSYRGADV FLLAFSLISK ASYENIHKKW LPELKHYPG IPIVLVGTKL DLRDDKQFLK DHPGAASITT AQGEELRKMI GAVRYLECSS KTQQNVKAVF DTAIRVALRP PKAKKKIKPL KTKRSRIC
Source	Yeast
Target Names	ARAC2
Protein Names	Recommended name: Rac-like GTP-binding protein ARAC2 Alternative name(s): GTPase protein ROP7
Expression Region	1-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.