



# Recombinant Arabidopsis thaliana Shaggy-related protein kinase eta (ASK7)

<b>Product Code</b>	CSB-MP655093DOA
<b>Abbreviation</b>	ASK7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q39011
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MADDKEMPAA VVDGHDQVTG HIISTTIGGK NGEPKQTISY MAERVVGTGS FGIVFQAKCL ETGETVAIKK VLQDRRYKNR ELQLMRVMDH PNVVCLKHCF FSTTSKDELF LNLVMEYVPE SLYRVLKHYS SANQRMPLVY VKLYMYQIFR GLAYIHNVAG VCHRDLKPQN LLVDPLTHQV KICDFGSAKQ LVKGEANISY ICSRFYRAPE LIFGATEYTT SIDIWSAGCV LAELLLGQPL FPGENAVDQL VEIIKVLGTP TREEIRCMNP HYTDFRFPQI KAHPPWHKIFH KRMPPEAIDF ASRLLQYSPS LRCTALEACA HPFFDELREP NARLPNGRPF PPLFNFKQEV AGSSPELVNK LIPDHIKRQL GLSFLNQSGT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ASK7
<b>Protein Names</b>	Recommended name: Shaggy-related protein kinase eta EC= 2.7.11.1 Alternative name(s): ASK-eta Protein BRASSINOSTEROID INSENSITIVE 2 Protein ULTRACURVATA 1
<b>Expression Region</b>	1-380
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.