



Recombinant Arabidopsis thaliana Dof zinc finger protein DOF3.4 (DOF3.4)

Product Code	CSB-MP655098DOA
Abbreviation	DOF3.4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q39088
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MPTSDSGEPR RIAMKPNGVT VPISDQQEQL PCPRCDSSNT KFCYYNNYNF SQPRHFCKAC RRYWTHGGTL RDVPVGGGTR KSAKRSRTCS NSSSSSVSGV VSNSNGVPLQ TTPVLFQSS ISNGVTHTVT ESDGKGSALS LCGSFTSTLL NHNAAATATH GSGSVIGIGG FGIGLGSGFD DVSFGLGRAM WPFSTVGTAT TTNVGSNGGH HAVPMPATWQ FEGLESNAGG GFVSGEYFAW PDLSITTPGN SLK
Source	Mammalian cell
Target Names	DOF3.4
Protein Names	Recommended name: Dof zinc finger protein DOF3.4 Short name= AtDOF3.4 Alternative name(s): OBF-binding protein 1
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.