



Recombinant Arabidopsis thaliana Serine--tRNA ligase (At5g27470)

Product Code	CSB-BP655107DOA
Abbreviation	At5g27470
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q39230
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MLDINLFREE KGNNPEIIRE SQRRRFASVE IVDEIIKLDK EWRQRQFEVD SFRKEFNKLN KQVAQLKIKK EDASEIIQQT EKNKQDSTAK EAEVREAYAA LKAKLEQVGN LVHDSVPVDK DEANNLVIKL WGEKRFSTPG LKLKNHVDLV ELLGIAADTKR GAEIAGARGF FLKGDGLMLN QALINFGLTF LKCRGFTGLQ PPFFMRKDVM AKCAQLAQFD EELYKVTGEG DDKYLIATAE QPLCAYHIDE WIHPTELPLR YAGYSSCFRK EAGSHGRDTL GIFRVHQFEK IEQFCITGPN ENASWEMLDE MMKNSDFYQ ALKLPYQIVS IVSGALNDAA AKKYDLEAWF PSSETFRELV SCSNCTDYQA RRLEIRYGQK KSNEQTKQYV HMLNSTLTAT ERTICCILEN YQREDGVDIP EVLQPFMGGE TFLPFKAKPV VADTKGKSK A
Source	Baculovirus
Target Names	At5g27470
Protein Names	Recommended name: Serine--tRNA ligase EC= 6.1.1.11 Alternative name(s): Seryl-tRNA synthetase Short name= SerRS Seryl-tRNA(Ser/Sec) synthetase
Expression Region	1-451
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.