



Recombinant *Synechococcus elongatus* Phosphate acyltransferase (plsX)

Product Code	CSB-BP655335FPY
Abbreviation	plsX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31N85
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MTRARIAVDA MGGDFAPEEI VKGALRAQEE LQADVILVGD PDRLRAICQD HAPHLQVRIE AEEAIAMED AAVSVRSRPR ASINVAMDV KAGEADAVIS AGHSGAVMAS ALLRLGRIRG IDRPAIGALL PTVIPGKPV LLDVGVANVDC KPRFLEQFAV MGSISRDVL GQANPRVGLV NIGEEDSKGN ELALASHQLL RNNPRICFVG NAEGRDVLSG QFDVVVCDGF VGNVLLKFAE AVGSVFLEII RDELPRGMRG KVGSTLLRRN LRRIKQRDLH AEHGGALLG VNGICIISHG SSKAPSIYSA IRLAVEAAEN RVIDHLHQIQ EPPSPAADLV TEPVVS
Source	Baculovirus
Target Names	plsX
Protein Names	Recommended name: Phosphate acyltransferase EC= 2.3.1.n2 Alternative name(s): Acyl-ACP phosphotransacylase Acyl-[acyl-carrier-protein]-phosphate acyltransferase Phosphate-acyl-ACP acyltransferase
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.