



Recombinant *Synechococcus elongatus* Methylenetetrahydrofolate--tRNA- (uracil-5)- methyltransferase TrmFO (trmFO)

Product Code	CSB-YP655337FPY
Abbreviation	trmFO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31NM4
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MAAISQPVIV IGAGLAGTEA AWQIAEAGVP VILYEMRPQR QSPAHHSESF AELVCSNSFG AMASDRAAGL LHEELRRLGS LVFSKASEHQ VPAGGALAVD RALFSEDLTR TVADHPLVEI RREELRSLPT EGIVVLCTGP LTSPDLAEDL QRFTGQDYCS FFDAASPIVT GESIDQAI AF RASRYDKGEA AYLNCP LN RD QYLAFREALV TAEQAELKDF EQESAKFFEG CLPIEELARR GEDTMRYGPL KPVGLFDARL GDWRDPENRS RRPYAI VQLR QEDRAGNLWN LVGFQTNLRW GEQKRIFQMI PGLSQAEFVR FGVMHRNTFV NAPQLLDASL QFRQRPTLLA AGQLIGTEGY SAAVAGGWLA GTNAARLALG RSPLVLPD TL VSGSLFRFIS SAEPKYFQPM PPNFGILPNL EQPPRNKKDR YAA YRDRALQ DLRDWQQQTHA IGNLSPSTGL ATTAIA
Source	Yeast
Target Names	trmFO
Protein Names	Recommended name: Methylenetetrahydrofolate--tRNA-(uracil-5)-methyltransferase TrmFO EC= 2.1.1.74 Alternative name(s): Folate-dependent tRNA (uracil-5)-methyltransferase Folate-dependent tRNA(M-5-U54)-methyltransferase
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.