



Recombinant Zea mays Uncharacterized protein ycf73 (ycf73-A)

Product Code	CSB-MP655693ZAX
Abbreviation	ycf73-A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q36907
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	≥85% (SDS-PAGE)
Sequence	MTKDETLLEV TLVSSVSIF LFGILLFMVL ISATRDFRER TSKLKVIMI WAGIVVITFA IAVRIYPIFI FLLKERIKPL VEALYDKLPW IWEVSLSRW DRLIDFLDRY LWACAQRIQT GIRKQKGEFV VTFSCRVKKR LYARAIEVGI HLSLLSNLFW ILKTTLAVGY RLL
Source	Mammalian cell
Target Names	ycf73-A;
Protein Names	Recommended name: Uncharacterized protein ycf73 Alternative name(s): ORF173
Expression Region	1-173
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.