



Recombinant *Synechococcus elongatus* 6,7-dimethyl-8-ribityllumazine synthase (ribH)

| | |
|--------------------------|---|
| Product Code | CSB-BP656017FPY |
| Abbreviation | ribH |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q31KZ5 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MAVFEGSFVN ASQFRLAVVI GRFNDLVCEK LLAGCQDCLK RHGVDPPDPQG TQVDYAWVPG SFEIPLVARQ LALTGRYDAI ICLGAVIRGQ TPHFDYVASE VAKGVAATSL QTGIPIAFGV LTVDNLQQAL ERAGIKSNLG WNYALSALEM ASLMHQIRST TGEVPRQTL P LAAAPSNFAG T |
| Source | Baculovirus |
| Target Names | ribH |
| Protein Names | Recommended name: 6,7-dimethyl-8-ribityllumazine synthase Short name= DMRL synthase Short name= Lumazine synthase EC= 2.5.1.9 Alternative name(s): Riboflavin synthase beta chain |
| Expression Region | 1-181 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |