



Recombinant *Synechococcus elongatus* Histidinol-phosphate aminotransferase (hisC)

Product Code	CSB-MP656034FPY
Abbreviation	hisC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31PF9
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	≥85% (SDS-PAGE)
Sequence	MLPFLRSELA RCQPYHPNPG GTGMAMDILD TNECPYDLPT DLKQTLADRY VEAIASNRYP DGSHTDLKAA IDYLSEQTA GQWQPGPEHV TVGNQSDELI RSILIATCLG GQGSVLVAEP TFSMYGIVAE TLGIPVVRIG RDPQTWEMDL AAAETAITQT EGTPVRLCFV VHPNSPTANP LTEAEKDWLR QVPPQILVVI DEAYFEFSGE TLLAELPQHP NWLITRTFSK ALRLAAHRVG YGIGDPQLIA ALEAIRLPYN LPSVAQLAAT LALEARSQLL SAIPRLITER DRLYRKLQVV SQLQVWPSAS NFLFLKTQSS SQTAAALAAQL KAQGTLVRHT ADGLRITIGS PAENERTLAH LQTAITQSLP ATV
Source	Mammalian cell
Target Names	hisC
Protein Names	Recommended name: Histidinol-phosphate aminotransferase EC= 2.6.1.9 Alternative name(s): Imidazole acetol-phosphate transaminase
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.