



Recombinant *Synechococcus elongatus* Phosphoribosylformylglycinamide synthase 1 (purQ)

Product Code	CSB-MP656037FPY
Abbreviation	purQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q31Q18
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MNVGIVVFPG SNCDRDVQWV TAGLLGQSTR MIWHEERDLS GLDLIVVPGG FSYGDYLRCG AIARFSPAMQ ATVAFAEAGG LVLGICNGFQ ILTEVGLLPG ALVRNRDLHF RCETTPLRVE RSDRPWSRTY QQGQILNLPI AHGEGRYHAD PATLALEAN GQVLFYLDN PNGSCNDIAG ITNVAGNVLG MMPHPERAAE AIAGSVDGLG LFAGLLEPVA A
Source	Mammalian cell
Target Names	purQ
Protein Names	Recommended name: Phosphoribosylformylglycinamide synthase 1 EC=6.3.5.3 Alternative name(s): Phosphoribosylformylglycinamide synthase I Short name= FGAM synthase I
Expression Region	1-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.