



Recombinant Danio rerio Serine/threonine-protein kinase pdik1l (pdik1l)

Product Code	CSB-MP656298DIL
Abbreviation	pdik1l
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q32PP3
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MVSSQAKYEL IQEVGRGSYG VVYEAVVRQT GARVAVEKIR CHSPENVELA LREFWALSSI QSQHPNVIHL EECVLQRDAL AQRMSHGSSS SLYLELVETS LKGEITFDPC CAYYMWFVMD FCDGGDMNAY LLSRKPSRKT NTSFMLQLGS ALAFLHRNQI IHRDLKPDNI LISQGRTPAG SPEPTLKVAD FGLSKVCSGS GLNPEEPASV NKCFLSTACG TDFYMAPEVW EGHYAKADI FALGVIIWAM VERITFVDVE TQKELLGSYV QQGEDIVPLG EALLENPKME LNIPARKKSM NASMKQLIRE MLSANPQERP DAFELELRLV RIACRELDWD T
Source	Mammalian cell
Target Names	pdik1l
Protein Names	Recommended name: Serine/threonine-protein kinase pdik1l EC= 2.7.11.1
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.