



Recombinant Human Transcription factor RFX4 (RFX4)

Product Code	CSB-BP656361HU
Abbreviation	RFX4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q33E94
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MHCGLLEEPD MDSTESWIER CLNESENKRY SSHTSLGNVS NDENEEKENN RASKPHSTPA TLQWLEENYE IAEGVCIPRS ALYMHYLDLC EKNDTQPVNA ASFGKIIRQQ FPQLTTRRLG TRGQSKYHYY GIAVKESSQY YDVMYSKKGA AWVSETGKKE VSKQTVAYSP RSKLGTLLPE FPNVKDLNLP ASLPEEKVST FIMMYRTHCQ RILDTVIRAN FDEVQSFLH FWQGMPPHML PVLGSSTVVN IVGVCDISLY KAISGVLMPV VLQALPDSLQ QVIRKFAKQL DEWLKVALHD LPENLRNIKF ELSRRFSQIL RRQTSLNHLC QASRTVIHSA DITFQMLEDW RNVDLNSITK QTLTYMEDSR DEHRKLITQL YQEFDHLLLE QSPIESYIEW LDTMVDRCVV KVAARKRQGS KKVAQQFLLM WSCFGTRVIR DMTLHSAPSF GSFHLIHLMF DDYVLYLLES LHCQERANEL MRAMKGEGST AEVREEIILT EAAAPTPSPV PSFSPAKSAT SVEVPPSPSS VSNPSPEYTG LSTTGAMQSY TWSLTYYTVT AAGSPAENSQ QLPCMRNTHV PSSSVTHRIP VYPHREEHGY TGSYNYGSYG NQHPPHMQSQ YPALPHDTAI SGPLHYAPYH RSSAQYPFNS PTSRMEPCLM SSTPRLHPTP VTPRWPEVPS ANTCYTSPSV HSARYGNSSD MYTPLTTRRN SEYEHMQHFP GFAYINGEAS TGWAK
Source	Baculovirus
Target Names	RFX4
Protein Names	Recommended name: Transcription factor RFX4 Alternative name(s): Regulatory factor X 4 Testis development protein NYD-SP10
Expression Region	1-735
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is a member of the regulatory factor X gene family, which encodes transcription factors that contain a highly-conserved winged helix DNA binding domain. This protein is structurally related to regulatory factors X1, X2, X3, and



X5. It has been shown to interact with itself as well as with regulatory factors X2 and X3, but it does not interact with regulatory factor X1. This protein may be a transcriptional repressor rather than a transcriptional activator. Three transcript variants encoding different isoforms have been described for this gene.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.