



Recombinant Enterobacteria phage T6 Uncharacterized 14.8 kDa protein in frd-Gp32 intergenic region (frd.2)

Product Code	CSB-MP656419EEA
Abbreviation	frd.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q38169
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T6 (Bacteriophage T6)
Purity	≥85% (SDS-PAGE)
Sequence	MEIGKKYELN PHRIKSFIDI SSSNANMVGI IQENGGWFEV KSISSLDGFD YVTEIICANG EIYNDGMDG DYFELSEEEF YCFREYKEPT SEEDEVKDKV SGVTKIHCIV DENNVDEIIE LLRKTFKK
Source	Mammalian cell
Target Names	frd.2
Protein Names	Recommended name: Uncharacterized 14.8 kDa protein in frd-Gp32 intergenic region
Expression Region	1-128
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.