



Recombinant Arabidopsis thaliana Heterodimeric geranylgeranyl pyrophosphate synthase small subunit, chloroplastic (GGR)

Product Code	CSB-EP656499DOA-B
Abbreviation	GGR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q39108
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	SSSSSAP GSLNFDLRTY WTTLITEINQ KLDEAIPVKH PAGIYEAMRY SVLAQQAKRA PPVMCVAACE LFGGDRLAAF PTACALEMVH AASLIHDDL CMDDDPVRRG KPSNHTVYGS GMAILAGDAL FPLAFQHIVS HTPPDLV PRATILRLITEIA RTVGSTGMAA GQYVDLEGGP FPLSFVQEKK FGAMGECSAV CGLLGGATE DELQSLRRYG RAVGMLYQVV DDITEDKKKS YDGGAEKGMM EMAEELKEKA KKELQVFDNK YGGGDTLVPL YTFVDYAAHR HFLPL
Source	E.coli
Target Names	GGR
Protein Names	Recommended name: Heterodimeric geranylgeranyl pyrophosphate synthase small subunit, chloroplastic
Expression Region	34-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.